

# Introducing a New Range Of Mobile Dental Units

## *Vet Pro ~ S and Vet Pro ~ B*

### These Units share our Goals for Quality & Reliability

|                     |  |
|---------------------|--|
| <b>CE</b> Marked    | The hall mark of Quality Assurance and conformity to European Safety Standards   |
| Proven Track Record | Manufacturer has been supplying Dental Units for use in human dentistry for many years - designs and longevity have been thoroughly road tested                                |
| Air Compressor      | Is from a major European compressor Manufacturer supported by a large English compressed air systems company   |
| Maintenance         | Technicians will install, instruct staff on maintenance procedures and annually Service these Dental Units & Compressors   |
| Confidence          | We and our collaborating manufacturers are confident that with our combined quality and maintenance regimes we can provide (subject to T & Cs) an outstanding value for money: |

### THREE YEARS PARTS WARRANTY - DENTAL UNITS & COMPRESSORS



## Vet Pro ~ B

Easy to operate Dental Unit - simple design manual selection of Polisher or Drill

Large operating table top area 270 x 370mm can easily accommodate an electric scaler

Water port for connection of Table Top Scaling Unit

Three way syringe provided

Ultra quiet semi housed compressor

### OPTIONS

- A. *Vet Pro ~ B* Mobile Unit with 3 way Syringe ultra quiet semi-housed compressor **£ 1,949**
- B. *Vet Pro ~ B* as A above with High Speed Turbine Handpiece, Air Motor, Slow speed handpieces - Straight & Propoxy for snap-on cups **£ 2,395**
- C. *Vet Pro ~ B* as B above and including Scaler incorporation kit to operate a table top scaler from single foot pedal - less clutter on floor **£ 2,489**

Please turnover page for details of the exciting *Vet Pro ~ S*