

Introducing a New Range Of Mobile Dental Units

Vet Pro ~ S

The elegant smooth lines of the Vet Pro ~ S is the hallmark of its overall characteristics

Surfaces are smooth and easy to clean ensuring good cross-infection control

Control knobs easily accessible yet concealed preventing ingress of bacteria and debris

Automatic switching between Drill & Polisher

Instrument Hanger Bar with Port for air scaler or table top scaler handpiece

Simple flushing system for all handpiece water channels to purge stagnant water to prevent legionella bacteria problems

Stainless Tray supplied for Instruments & Table Top Scaler

Upgrades for Fibre Optics and other bespoke options



OPTIONS

- A. **Vet Pro ~ S Mobile Unit with 3 way Syringe ultra quiet semi-housed compressor** **£ 2,399**
- B. **Vet Pro ~ S as A above with High Speed Turbine Handpiece, Air Motor, Slow speed handpieces - Straight & Prophy for snap-on cups** **£ 2,845**
- C. **Vet Pro ~ S as B above and including Scaler incorporation kit to operate a table top scaler from single foot pedal - less clutter on floor to trip over** **£ 2,930**

(All prices are exclusive of VAT)