

MAGW100

Wood's light

Technology that's clear to see

The **DARAY** MAGW100 Wood's light is a must-have diagnostic tool for analysing skin conditions.


The MAGW100 allows clinicians to correctly inspect and diagnose a variety of conditions. Lightweight and easy to use, the MAGW100 emits an ultraviolet light with the remarkable ability to note skin colour changes and variations in skin tone.



- Emits an ultraviolet light allowing clinicians to correctly inspect and diagnose a variety of conditions.
- Remarkable ability to note skin colour changes and variations in skin tones.
- Used in a variety of applications including dermatology, forensic science, veterinary medicine and childcare.
- Lightweight and easy to use

Specifications

4 x 16W flourescent lamps	Protective cover plate
~5,000 hrs life per lamp	230V mains operated
5x magnification	Dimensions 330x195x70mm



Percent irradiance



ABOUT DARAY

DARAY have been established in the medical, dental and veterinary markets since 1968 building exceptional and durable products.

We shall continue to innovate; aiming to excel in our endeavours to deliver reliable solutions.



DARAY
HEALTHCARE PRODUCTS

www.daray.com

Tel: 0800 804 8384 Email: info@daray.co.uk