



AGROSCAN A14

DISPLAY

5.5 High resolution LCD
256 grey levels

WEIGHT

1.5 kg (battery included)

SIZE

9x23x10.5 cms

PROBE

- Mechanical sector of 3.5/5 MHz
- Special crystal angle to be used rectally on large animals

IMAGE MODE

B-mode
M-mode
Exploration : 3 depth
(8-11-15 cms)

BATTERY

- Inbuilt battery – memory effect
- Working time : more than 3 hours and 30 minutes

MEASUREMENT

- Distance
- Length
- Volume
- area & circumference
- ellipse
- back fat on small ruminats and swine

CONNECTION


- Standard video output (NTSC signal)

STANDARD ACCESSORIES

- Transportation case
- Working carrier
- Carrying straps
- Charger and cables
- instruction manual
- 250 ml gel

OPTIONS


- Sun protector
- Video goggles
- multimedia player
- ultrasound gel




Pregnancy 30 days on cow




Pregnancy 85 days on cow



Empty uterus on goat



Pregnancy 75 days on goat



Pseudo pregnancy on goat

