

# MICROMIX 80/40



## A CONTROL PANEL

POWER SUPPLY:	230V 50Hz
COUNTER POWER:	1.6KW
kV ADJUSTMENT:	80 kV
mA ADJUSTMENT :	40mA
sec. ADJUSTMENT:	from 0.01 to 1.99 sec

## B X RAY TUBE-HEAD

VOLTAGE:	80kV
CURRENT:	40mA
LEAKAGE RADIATION:	less than 1mG/hour according to IEC601-1-3

## C X RAY TUBE

FOCUS:	1.8mm
MAX VALUE OF KV:	90 kV
MAX VALUE OF mA:	120mA
ANOD MATERIAL:	Tungsten
ANOD ANGLE:	19°

## D SAFETY DEVICES

VISUAL INDICATOR TO SEE X RAY UNIT ALREADY PERFORMED

## E OPTIONAL

FOOT-SWITCH WITH 6m CABLE

