

V60

Handheld pulse oximeter
Technology that doesn't miss a beat



Computer not supplied

The **DARAY** V60 hand-held pulse oximeter is lightweight, portable and easy to use.

- Compact, light and easy to carry
- Clear, bright LED display
- Powered by 2 x AA batteries
- Data storage up to 24 hours
- Weight 136g (with batteries)
- Transfer data to PC for analysis, saving or printing
- Pulse waveform display
- Pulse rate: accuracy $\pm 2\%$ or $\pm 2\text{BPM}$
- Battery charge indicator

ABOUT DARAY

DARAY have been established in the medical, dental and veterinary markets since 1968 building exceptional and durable products.

We shall continue to innovate; aiming to excel in our endeavours to deliver reliable solutions.



FM98457

DARAY[®]
HEALTHCARE PRODUCTS

www.daray.com

Tel: 0800 804 8384 Email: info@daray.co.uk